



Routine Maintenance for Atomic Absorption Spectrophotometers

Application Note

Atomic Absorption

Author

Margaret A. Cunliffe

Introduction

Instruments in good operating condition are a necessity in any analytical laboratory. This level of integrity can be achieved by a regular maintenance schedule with minimal work. The four main areas of such a program for atomic absorption spectrophotometers include:

- General instrument maintenance
- Gas supply maintenance
- Flame component maintenance
- Furnace component maintenance

The benefits of routine maintenance include:

- Increased instrument lifetime
- Reduced downtime
- Overall improvement in instrument performance; giving the operator greater confidence in the validity of his analytical results



Agilent Technologies

General Instrument Maintenance

Dust and condensed vapors can accumulate on the instrument case, and corrosive liquids can be spilled on the instrument. To minimize damage, wipe off the instrument with a damp, soft cloth using water or a mild detergent solution. **DO NOT USE ORGANIC SOLVENTS.** The sample compartment windows and the lamp windows can accumulate dust or fingerprints. In such cases, clean the windows with a soft tissue moistened with a methanol or ethanol and water solution. If the windows are not clean, the operator will observe noisy lamp signals and non-reproducible analytical results.

The remaining optical components are sealed, but they should not be exposed to corrosive vapors or a dusty atmosphere. In laboratories where high concentrations of dust or vapors are unavoidable, schedule a yearly check by a service engineer to maintain the efficiency of optical light transmission in the instrument. There is no need for an operator to clean the sealed optical components.

Gas Supply Maintenance

Three gases are suitable for flame M. Air and nitrous oxide are used as combustion support gases (oxidants). Acetylene is used as the fuel gas. Each gas is supplied to the instrument through piped supply systems and rubber hoses. Copper or copper alloy tubing may be used for the oxidant gases. Acetylene should only be supplied through stainless steel or black iron pipe. Check connections regularly between the supply and instrument for leaks, especially when tanks are changed using a soap solution or commercial leak detector. Check the rubber hoses connected to the instrument for fraying and cracking. In addition, each time a tank is changed, check the regulators and valves for proper operation.

Because potentially toxic gases are used or produced in the flame, it is necessary to use a suitable exhaust system with a minimum capacity of 6 m³/min (200 cfm). A simple smoke test will indicate if it is functioning properly.

Compressed Air Supply

Air may be supplied to the instrument from cylinders, a house air system, or small compressor. Cylinders are the most expensive source of air, particularly where large amounts are consumed and cylinders must be changed frequently. If compressed air from an in-house supply is used, a filter/regulator assembly must be installed in the input line to the instrument. An acceptable "Air Service Unit" (Part No. 01 102093 00) may be ordered from any Agilent sales office.

Whatever source is used, the supply must be continuous and have a delivery pressure of 420 kPa (60 psi). The air must be clean, dry and oil free. Approximately 50% of all gas unit failures are caused by moisture or other impurities in the air supply.

Excessive noise in the readout has also been attributed to contaminated air. An air filter assembly is therefore an essential component of the atomic absorption spectrophotometer, and its inclusion in the air supply installation is mandatory. Weekly, check the air filter for particle and moisture accumulation. When necessary, dismantle the air filter assembly and clean the filter element, bowl, and drain valve components. Use the following procedure for dismantling and cleaning the air filters supplied with the instrument.

1. Shut off the air supply and allow the system pressure to bleed off.
2. Unscrew the filter bowl, complete with automatic drain valve.
3. Unscrew the retaining ring and push the drain valve back into the bowl.
4. Unscrew the baffle carefully, and remove the filter and filter shield.
5. Clean the filter bowl, drain valve components, baffle, and filter shield by washing in a solution of soap and water. **DO NOT USE ORGANIC SOLVENTS AS THEY WILL DESTROY THE BOWL AND VALVE COMPONENTS.** Rinse thoroughly in fresh water.
6. Clean the filter element by washing in ethyl alcohol or similar solvent.
7. Ensure that all components are properly dried before reassembly.

Nitrous Oxide Supply

The nitrous oxide used for atomic absorption spectrophotometry must be oil free. If a heated regulator is not used, loss of regulation can occur due to the expansion cooling effect encountered when nitrous oxide is drawn from a cylinder. This can lead to erratic results and create a potential flashback situation with manual gas control units: An acceptable heated regulator may be ordered from any Agilent sales office. The consumption rate is dependent on the application, but is usually 10–20 liters per minute.

Acetylene Supply

Acetylene is the only combustible gas which is normally used in MS. The gas must be supplied packed in acetone. Some companies supply acetylene packed in proprietary solvents, but unfortunately the disadvantages outweigh the advantages. The major disadvantage is that the solvent may be carried over into the instrument and corrode the internal tubing, causing a potential explosion hazard. Ensure that the acetylene is at least 99.6% pure "M Grade" and packed in acetone.

The delivery pressure must be regulated and never exceed 105 kPa (15 psi). Check the instrument operation manual for the correct delivery pressure for the particular instrument being used. In addition, check the acetylene cylinder pressure daily, and maintain in excess of 700 kPa (100 psi) to prevent acetone from entering the gas line and degrading analytical results or causing damage to the instrument.

Flame Component Maintenance

The flame component section of the instrument can be divided into three areas; the nebulizer, spray chamber and burner. Each requires routine maintenance to assure optimum performance.

Nebulizer

The nebulizer area of the flame component consists of the capillary tubing and the nebulizer body. Always ensure that the plastic capillary tubing used for aspirating solutions is correctly fitted to the nebulizer capillary. Any leakage of air, tight bends, or kinks will cause unsteady, non-reproducible readings.

At times the plastic capillary tubing can become clogged and it will be necessary to cut off the clogged section or fit a new piece of capillary tubing (about 15 cm long). In any event, make sure the plastic capillary tubing fits tightly on the nebulizer capillary. The nebulizer capillary can also become clogged. If this occurs, proceed as follows:

1. TURN THE FLAME OFF.
2. Remove the plastic capillary tubing from the nebulizer.
3. Remove the nebulizer from the bung.
4. Dismantle the nebulizer as described in the instrument operation manual or the instruction manual supplied with the nebulizer.
5. Place the nebulizer in an ultrasonic cleaner containing 0.5% liquid soap solution such as Triton X-100 for 5 to 10 minutes. If the ultrasonic bath fails to clear the block-

age, pass a burr-free nebulizer wire CAREFULLY through the nebulizer and then repeat the ultrasonic cleaning procedure.

6. Re-assemble the nebulizer in accordance with the instructions.
7. Install the cleaned nebulizer.

Replace the plastic capillary tubing.

If blockages are allowed to build up and are not removed, the analytical signal will steadily drop until no absorbance is observed.

8. Check the nebulizer body, capillary, and venturi occasionally for corrosion. Nebulizer problems can be minimized by taking care to always aspirate 50–500 mL of distilled water at the end of each working day.

Spray Chamber

As the sample leaves the nebulizer it strikes the glass bead and breaks into an aerosol of fine droplets. The efficiency of the glass bead can be degraded by surface cracks, pitting and the accumulation of solid material. The reduction in bead efficiency can cause lower absorbance readings and noisy signals. When removing the nebulizer for inspection, always check the glass bead. Look for pitting, cracks, breakage, ensure that the adjusting mechanism operates properly and that the bead is correctly positioned over the nebulizer outlet (venturi).

While the nebulizer and glass bead are removed from the instrument for inspection, the spray chamber and liquid trap should be removed, dismantled, and cleaned. Discard the liquid in the liquid trap and wash both the spray chamber and liquid trap thoroughly with laboratory detergent and warm water. Rinse completely with distilled water and dry all components. Refill the liquid trap and reassemble the spray chamber, checking for any distortion of O-rings or blockages in the gas inlets. Reconnect the drain hose. If a bottle or jug is used to collect the waste solutions, check that the hose is not below the level of the waste. If the hose is below that level, absorbance readings will steadily decrease with occasional abrupt increases as intermittent drainage of the spray chamber occurs. Therefore, it is necessary to daily check the level of the waste and to dispose of it frequently. This is imperative when using organic solvents because of the potential hazards introduced by flammable liquids. Only wide necked, plastic containers can safely be used to collect the waste solutions.

Burner

The final area of concern in the flame component is the burner. During aspiration of certain solutions, carbon and/or salt deposits can build up on the burner causing changes in

the fuel/oxidant ratio and flame profile, potential clipping of the optical beam, and degradation of the analytical signal. To minimize the accumulation of salts, a dilute solution of acid (HNO_3) may be aspirated between samples. However, if salts continue to build up, turn off the flame and use the brass cleaning strip supplied with the instrument. Insert the strip in the burner slot and move it back and forth through the slot. This should dislodge any particles which will then be carried away once the flame is lit and water aspirated.

DO NOT USE SHARP OBJECTS such as razors to clean the burner as they can nick the slot and form areas where salt and carbon can accumulate at an accelerated rate.

If this type of cleaning is inadequate, remove the burner, invert, and soak it in warm soapy water. A scrub brush will facilitate cleaning. Soaking may also be done in dilute acid (0.5% HNO_3). Ultrasonic cleaners containing dilute non-ionic detergent only are another alternative for cleaning. After cleaning, thoroughly rinse the burner with distilled water and dry before installing in the instrument. **NEVER DISASSEMBLE THE BURNER FOR CLEANING. IMPROPERLY RE-ASSEMBLED BURNERS WILL LEAK COMBUSTIBLE GAS MIXTURES, POTENTIALLY CAUSING EXPLOSIONS.**

Each day after all analyses are completed, 50–100 mL of distilled water should be aspirated to clean the nebulizer, spray chamber, and burner. This is even more important after aspirating solutions containing high concentrations of Cu, Ag, and Hg, since these elements can form explosive acetylides. The entire burner/nebulizer assembly should be disassembled and thoroughly cleaned after analyzing these types of solutions. The burner should be removed weekly, scrubbed with a laboratory detergent, and rinsed with distilled water.

Furnace Component Maintenance

The graphite furnace accessory maintenance can be divided into three major areas; the gas and water supplies, the workhead, and the autosampler. Each plays an important role in obtaining valid analytical results. The following general maintenance program refers to the GTA-95.

Gas and Water Supplies

Normally the gases used in FAAS are inert gases such as N_2 and Ar. Either one may be used, but must be clean, dry, and of high purity. The regulated pressure should be 100–340 kPa (15–50 psi). At times the incorporation of air may be useful to fully ash a sample. However, air should not be used at ash temperatures higher than 500 °C because of the accelerated rate of graphite component deterioration at elevated temperatures.

The water supply, used to cool the furnace, may be supplied either from a laboratory tap or a cooling-recirculating pump. If a recirculating pump is used the water must be kept below 40 °C. The water used must be clean and free of corrosive contamination. The flow should be 1.5–2 liters/minute. Maximum permissible pressure is 200 kPa (30 psi).

Workhead

The workhead is a closed assembly with quartz windows on either end. Before starting an analysis, check the windows for dust or fingerprints. If needed, clean both sides of the quartz windows with a soft tissue moistened with an alcohol/water solution. Never use coarse cloths or abrasive cleaning agents. While the windows are removed, inspect the gas inlets on the window mountings. If the graphite components have deteriorated extensively, graphite particulates may have dropped into the gas inlets, blocking the proper flow of gas. This will cause further graphite deterioration at an accelerated rate and lead to poor analytical performance. To clean, carefully blow out the particulates with a supply of air. Inspect the inside of the window mountings and clean off any sample residue which may have deposited over time.

In the center of the workhead are the graphite components. At frequent, regular intervals, remove the graphite tube atomizer and inspect the inside of the graphite shield. Ensure that the bore and the injector hole area are free of loose carbon or sample residue. Check the electrodes on either end of the graphite shield for proper tapering. If the tapering is worn or burnt, the electrodes will not make the correct contact with the graphite tubing, causing fluctuations in applied power resulting in irreproducibility. The electrodes also have a series of gas inlets which must be free of loose carbon or sample residue.

Above the graphite shield is the titanium chimney. Injected sample or sample residue from the ash/atomize cycles may deposit in this area. A cotton swab soaked with alcohol can be used to clean both the inside and outside of the chimney. Alternatively, the titanium chimney may be soaked in dilute acid to remove deposits.

Autosampler

The components of the autosampler requiring routine maintenance are the rinse bottle, syringe, and capillary tubing, the proper care of which will minimize contamination and improve reproducibility of analytical results.

Regularly remove the rinse bottle for cleaning. This involves soaking the bottle in 20% HNO_3 followed by rinsing with distilled-deionized water. Refill the bottle with a solution of 0.01–0.05% HNO_3 in distilled-deionized water. The solution

may also include 0.005% v/v Triton X-100 R. The Triton helps maintain the sample capillary in clean condition and assists in obtaining good precision.

At times, graphite particulates may accumulate on the capillary tip and should be carefully removed with a tissue. If these particulates are not removed, the dispensing characteristics of the capillary may change. Contamination of the capillary may become a problem when using some matrix modifiers. In such cases, direct the capillary to a vial containing 20% HNO₃, draw up 70 µL, and stop the autosampler while the capillary is in the vial. After a period of a few minutes, the autosampler RESET should be utilized to rinse out the acid solution. This will clean the internal and external areas of the capillary. Similarly, organic residues can be removed by directing the capillary to a vial of acetone and repeating the above procedure. The PTFE capillary should be treated carefully during cleaning and operation. If bends or kinks appear, it can take time to reshape, and while doing so the repeatability of injection may be degraded. If the capillary tip is damaged, the damaged portion should be cut off at a 90° angle with a sharp scalpel or razor blade.

The final area of the autosampler maintenance schedule is the syringe. Daily, check for bubbles in both the capillary and syringe. Any bubbles in the system can cause dispensing errors and lead to erroneous results. Follow the instructions in the operating manual to free the system of bubbles. If the bubbles continue to cling to the syringe, it may need cleaning. The syringe can be washed with a mild detergent solution and thoroughly rinsed with deionized water. Ensure that contamination is not introduced through the syringe. Be particularly careful not to bend the plunger while washing the syringe.

Conclusion

Attached is a routine maintenance schedule for atomic absorption spectrophotometers (Figure 1). By adhering to this program, the overall integrity of the atomic absorption spectrophotometer can be maintained and the laboratory analyst will reap the benefits of increased instrument lifetime, reduced downtime, and gain greater confidence in the analytical results.

| Maintenance Schedule (Flame AA) | |
|---|------------------|
| Daily | Completed |
| 1. Check Gas | |
| 2. Check Exhaust system with smoke test | |
| 3. Empty the drain receptacle | |
| 4. Clean lamp and sample compartment windows | |
| 5. Rinse spray chamber with 50-100 mL of distilled water | |
| Weekly | |
| 1. Disassemble spray chamber | |
| (a) Check glassbead | |
| (b) Check nebulizer components | |
| (c) Wash the spray chamber and liquid trap | |
| (d) Scrub the burner | |
| (e) Change the liquid in the liquid trap | |
| (f) Check the O-rings | |
| 2. Check air filter assembly | |
| 3. Wipe off instrument | |
| 4. At Time of Gas Tank Change | |
| 5. Check for leaks | |
| 6. Check for operation of the regulators | |
| 7. Check for operation of the shut off valves | |
| 8. Check the gas supply hoses | |
| Yearly | |
| 1. Schedule an Agilent service engineer to perform Preventive Maintenance | |

Figure 1. Routine maintenance schedule for atomic absorption spectrophotometers.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 1984
 Printed in the USA
 November 1, 2010
 AA039

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem



Agilent Technologies